

# ISD-15-BTS2048-VL

<https://www.gigahertz-optik.com/en-us/product/isd-15-bts2048-vl/>


**Product tags:** VIS



# Description

## The First Class BTS2048 Spectroradiometer Series

The high-quality [BTS2048 Series](#) CCD based spectroradiometers are internationally recognized as a high-end products. It is one of the most compact spectroradiometers on the market which enables direct system integration in many applications without the need for expensive, and potentially measurement degrading, light guides. Among its characteristic features is its diffuser window which permits direct mounting onto integrating spheres (e.g. the ISD-15) for measurement of luminous flux. More detailed information about the [BTS2048-VL](#) can be found in the respective data sheets. The unit is also available in the [BTS2048-VL-TEC](#) thermoelectric cooled version.




## Use of the BTS2048 Series for front-end and back-end LED binning

The BTS2048 is ideally suited for use in industrial LED binning. In particular, it can be used to perform time-synchronized, pulsed measurements of LEDs as stipulated by CIE S025 or DIN 5032-9. This is necessary when one wants to account for thermal effects in measurements (see also technical article on [SSL/LED testing](#)). The instantaneous, electronic zero setting of all pixels (electronic shutter) of the CCD sensor in conjunction with its trigger interface, results in perfect synchronization between the BTS2048 measurement and the LED current supply. The device also has an exceptionally wide dynamic range with a minimum integration time of just 2  $\mu$ s. This, is 1000 times shorter than conventional devices (typically milliseconds), hereby makes it possible to achieve an OD3 factor. Additionally, the internal OD1 and OD2 filters increase the dynamic range up to 9 orders of magnitude and the fast Ethernet interface provides rapid data transfer. For users to integrate the device into their system, the BTS2048 Series is supplied with the powerful [S-SDK-BTS2048](#) software development kit.

- 1) BTS2048-VL
- 2) BiTec sensor with Si photodiode, CCD spectroradiometer
- 3) Filter wheel with OD1, OD2 and shutter
- 4) USB 2.0 Interface
- 5) High speed ethernet Interface
- 6) Trigger In/Out
- 7) Diffuser in sphere surface level
- 8) Measurement port with protection dome window and knife-edge design
- 11) Auxiliary lamp
- 12) Baffle

\*One of its unique features is the from Gigahertz-Optik developed innovative [BiTec sensor](#) that consists of a  $V(\lambda)$  filtered Si photodiode and a spectroradiometer unit. This makes it extremely linear, stable, and fast and is therefore a guarantee for higher measurement accuracy which is not accompanied by any disadvantages. Both sensors can be used independently and the mutual correction of the sensors is advantageous for accuracy, speed and versatility (see article on [BTS technology](#)).




## Enhancement of the BTS2048 to a 2Pi luminous intensity spectroradiometer

The BTS2048 can be combined directly with integrating spheres in order to measure luminous intensity. With its compact design, light weight, and no need of a connecting light guide, the integrating sphere and spectroradiometer form a monolithic module that can be conveniently integrated in complete test systems for front-end and back-end LED binning. It is supplied as a fully calibrated module. The module can be recalibrated as a single entity. Eliminating the need for connecting light guides enables convenient recalibration with low uncertainty. The [BN-LHSE-2P-20](#) 2pi calibration lamp can be used to perform the recalibration directly in the test system using the supplied user software.

## S-BTS2048 application software

The [Integrating Sphere Construction Kit](#) from Gigahertz Optik GmbH allow for customization and hence the best possible adaptation of the integrating sphere to the specific application.



## The ISD-15-V01 integrating sphere for industrial LED binning

The ISD-15-V01 has many properties that make it perfect for luminous flux measurements in combination with the BTS2048 spectroradiometer. It has a 150 mm (6 inch) diameter and offers a good compromise in terms of 2Pi acceptance and light throughput for high sensitivity. The 35 mm diameter measurement port has a knife edge for unrestricted light entry. A transparent dome behind the measurement port is integrated to protect the coating of the sphere from contamination in industrial environments. The sphere features an auxiliary lamp to compensate for the self-absorption effects. Illumination with the auxiliary lamp is done through a diffuser window.

## S-SDK-BTS2048 Software Development Kit

## Calibration

One essential quality feature of photometric devices is their precise and traceable calibration. The ISD-15-V01 -BTS2048-VL 2pi luminous flux spectroradiometer is calibrated at Gigahertz-Optik's DAKS accredited (D-K-15047-01-00) [ISO/IEC 17025 calibration laboratory](#) for the *spectral responsivity* and *spectral irradiance*. The calibration is performed with a [BN-LHSE-2P-20](#) calibration lamp that has 2pi radiation characteristics in the integrating sphere. Every device comes with its respective calibration certificate.

## Software for system integration

Gigahertz-Optik offers the S-SDK-BTS2048 software development kit for integration of the

## Specifications




### General


|                      |                                                                                                                                                                                          |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Short description    | 2pi luminous flux spectroradiometer as an integration module for LED test systems for front-end and back-end LEDs based on the luminous flux, spectrum, color, and Color rendering index |
| Main features        | Integrating sphere with a 15 cm diameter and a 35 mm measurement port that is protected against contamination. High-quality CCD sensor spectroradiometer                                 |
| Measurement range    | Luminous flux: Integral 0.03 mlm to 100 klm, Spectral for typical white LEDs: 2 lm to 16000 klm<br>Spectral radiant flux: 1E-7 W/nm to 7E2 W/nm<br>Spectral range: 350 nm to 1050 nm     |
| Typical applications | 2pi luminous flux spectroradiometer for integration in LED test systems. In-line quality assurance of LED circuit boards, LED matrices, LED spot lamps                                   |
| Calibration          | Factory calibration. Traceable to international calibration standards                                                                                                                    |

## Downloads

| Type                    | Description                   | File-Type | Download                                                                                                                                                                          |
|-------------------------|-------------------------------|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| BTS2048-Series Brochure | Not Just Another Spectrometer | pdf       | <a href="https://www.gigahertz-optik.com/assets/BTS2048_broschuere_DI_NA4_hoch_V2_2022.pdf">https://www.gigahertz-optik.com/assets/BTS2048_broschuere_DI_NA4_hoch_V2_2022.pdf</a> |

## Configurable with

| Product Name   | Product Image                                                                       | Description                                                                                        | Go to product                                                                                                                             |
|----------------|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| BTS2048-VL     |  | Versatile superb Speed and high Quality LED Spectroradiometer                                      | <a href="https://www.gigahertz-optik.com/en-us/product/bts2048-vl/">https://www.gigahertz-optik.com/en-us/product/bts2048-vl/</a>         |
| BTS2048-VL-TEC |  | Versatile Temperature Controlled High Speed and High Quality LED Spectroradiometer                 | <a href="https://www.gigahertz-optik.com/en-us/product/bts2048-vl-tec/">https://www.gigahertz-optik.com/en-us/product/bts2048-vl-tec/</a> |
| BTS2048-IR     |  | Compact IR Spectroradiometer Fulfilling all the Requirements of a High-End Array Spectroradiometer | <a href="https://www.gigahertz-optik.com/en-us/product/bts2048-ir/">https://www.gigahertz-optik.com/en-us/product/bts2048-ir/</a>         |

| Product Name   | Product Image                                                                     | Description                                                                                                                        | Go to product                                                                                                                             |
|----------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| BTS2048 Series |  | Compact spectroradiometers with excellent optical performance and BiTec technology for precise measurements for lab and field use. | <a href="https://www.gigahertz-optik.com/en-us/product/bts2048-series/">https://www.gigahertz-optik.com/en-us/product/bts2048-series/</a> |

## Purchasing information

| Article-Nr         | Modell                 | Description                                                                                                                                                                                                                                    |
|--------------------|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b>     |                        |                                                                                                                                                                                                                                                |
| 15306340           | ISD-15-V01             | Integrating sphere, entrance port with dome window, auxiliary lamp with baffle, detector port with UMPF-1.0-HL and baffle, port caps.                                                                                                          |
| 15298281           | BTS2048-VL             | Measuring device, hard cover box, users guide, S-BTS2048 software, calibration certificate.                                                                                                                                                    |
| 15298687           | BTS2048-VL-TEC         | Measuring device, hard cover box, users guide, S-BTS2048 software, calibration certificate.                                                                                                                                                    |
| <b>Calibration</b> |                        |                                                                                                                                                                                                                                                |
| 15300771           | K-BTS2048VL-Phi2-S-V01 | Calibration of the spectral radiant flux in W/nm from 350 nm to 1050 nm and of the luminous flux in lm of the BTS2048-VL with optional integrating sphere. Use of a calibration standard with 2Pi light distribution. Calibration certificate. |
| <b>Software</b>    |                        |                                                                                                                                                                                                                                                |
| 15298470           | S-SDK-BTS2048          | Software development kit, software CD with users guide.                                                                                                                                                                                        |

## Contact, Calibration, Service & Support

We are known worldwide for excellent technical consulting and after sales support. Contact us to find together the best solution for you. Our services:

- Technical Consulting & Sales
- After-Sales Support
- Calibrations & Re-Calibrations ([ISO/IEC 17025 Calibration Services](#), [factory calibration](#), [Calibration of Third-Party Products](#))
- Repairs & Updates
- OEM & Feasibility Consulting of Customized Solutions

[Send us your inquiry](#) or contact us by phone or e-mail. We would welcome your feedback too or review us on [Google](#).

### Gigahertz Optik GmbH (Headquarter)

Tel.: +49 (0)8193-93700-0  
Fax: +49 (0)8193-93700-50  
[info@gigahertz-optik.de](mailto:info@gigahertz-optik.de)

An der Kaelberweide 12  
82299 Tuerkenfeld, Germany

### Gigahertz-Optik, Inc. (US office)

Phone: +1-978-462-1818  
[info-us@gigahertz-optik.com](mailto:info-us@gigahertz-optik.com)

Boston North Technology Park  
Bldg B - Ste 205  
Amesbury, MA 01913 USA